



# Waterloo

## HEALTHCARE

PATIENT BIN CART SYSTEMS



800:833:4419

direct 602:414:3691

[waterloohealthcare.com](http://waterloohealthcare.com)

WATERLOO HEALTHCARE

3730 East Southern Avenue  
Phoenix, Arizona 85040-3958

# MEDTRX & HALFTRX Patient Cassette Cart Systems

Select one of these configurations or build your own custom cart, see page 9. All carts shown with optional narcotics lock (NL-MT) and drawer trim (see page 9).

Antimicrobial properties built-in to protect portions of the product.



Antimicrobial Surfaces



MTK1-4B411/6TTTT



MTK1-4B411/TFTF6



MTK1-4B771/TFTFT



MTK1-4B771/33FTF6



MTK1-4B771/33TTTT




HTK1-6BTTTT



HTK1-33TTT6




## Medtrx/Halftrx Cassettes and Bins



### CASS-3MTC

Two-tier cassette includes (10) 3" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-3G)

16.5"W x 13.5"D x 6"H



### BIN-3G

Patient bins for MEDTRX medication cassettes (10/ctn).

2 7/8"W x 10 3/4"D x 2 3/16"H



### LBL-1B

Bin labels for MEDTRX bins. 1000 labels. Size: 2.5"W x 0.75"H



### CASS-5MTC

Two-tier cassette includes (6) 5" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-5G)


16.5"W x 13.5"D x 6"H



### BIN-5G

Patient bins for MEDTRX medication cassettes (6/ctn).


4 13/16"W x 10 3/4"D x 2 3/16"H



### CASS-7MTC

Two-tier cassette includes (4) 7" wide bins for MEDTRX series carts, exchange programs or individual use. Additional bins sold separately (BIN-7G)

16.5"W x 13.5"D x 6"H



### BIN-7G

Patient bins for MEDTRX medication cassettes (4/ctn).

7 1/4"W x 10 3/4"D x 2 3/16"H


## DMC Cassettes and Bins



### CASS-2S


Two-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-2S110

21"W x 10"D x 7.75"H



### MBG-3B

3" wide patient bins for DMC medication cart (24/ctn). 3"W x 9"D x 2.75"H. Dividers sold separately.



### BDM-2

Left to right bin dividers that fit 4.5" DMC bins. Fits part number MBG-4B.


Bins sold separately.



### CASS-3S


Three-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-3S110

21"W x 10"D x 10.75"H



### MBG-4B

4" wide patient bins for DMC medication cart (16/ctn). 4.5"W x 9"D x 2.75"H. Dividers sold separately.



### BDS-1


Left to right bin dividers that fit 3" DMC bins. Fits part number MBG-3B. Bins sold separately.



### CASS-4S


Four-tier lockable cassette for dual-sided medication carts, exchange programs or individual use. Bins sold separately. Multiple orders can be "keyed alike" just use part number CASS-4S110

21"W x 10"D x 13.75"H



### MBG-6B

6" wide patient bins for DMC medication cart (12/ctn). 6"W x 9"D x 2.75"H. Dividers sold separately.



### BDW-6


Left to right bin dividers that fit 6" DMC bins. Fits part number MBG-6B. Bins sold separately.



### PBL-1


Bin labels for plastic bins in DMC carts (MBG-3B, 4B, 6B). 1 roll of 1000 labels for individual patient bins.

## Medtrx/Halftrx Drawer Dividers




### DIV-HTS

Divides half-size 4" - 7" drawers from front to back. Includes dividers and brackets (3.5"H).




### DIV-HTL

Divides half-size 9" or larger drawers from front to back. Includes dividers and brackets (8.5"H).



### DIV-MTS

Divides full-size 4" - 7" drawers from front to back. Includes dividers and brackets (3.5"H).



### DIV-MTL

Divides full-size 9" or larger drawers from front to back. Includes dividers and brackets (8.5"H).



## Cart Accessories



**PC1000**  
1000 Pill crushing replacement pouches. Tablets are crushed in a durable latex-free pouch.




**BSK-1**  
Small (6 quart) plastic waste basket. Requires mounting rail (SR-1).  
13"W x 6.3"D x 10"H



**SB-2**  
Non-locking sharps container mounting bracket. Mounts to side of all, Medtrix and DMC. Designed for 5 quart container.



**WL-3**  
Screws into wall and allows you to secure your medication cart to the wall.



**BSK-2**  
Large plastic waste basket. Requires mounting rail (SR-1). 13"W x 8.5"D x 15"H



**SB-2L**  
Sharps bracket includes key lock for added safety & security. Fits 5 quart sharps container.



**CDA-2**  
Adjustable pop-up style cup dispenser. Holds 2 to 5 oz. cups. For use with all of our carts.



**SR-1**  
A side rail for hanging accessories on sides of carts.



**LH-1**  
Label holders, self-adhesive. 1"H x 3"W, package of 100. Holders are loaded from the side.



**EWR-1**  
Wire trash container frame with lid and 10 bags. Replacement bags available (EWR-BAG)  
12.5"W x 7"D x 19.5"H



**SB-1**  
Sharps container mounting bracket with adjustable velcro strap. Mounts to side of all Medtrix and DMC.

## Cart Accessories



**OS-1**  
Hospital grade outlet strip. 15 foot long cord with 6 outlets.



**IV-2**  
Chrome IV pole and mounting bracket. Mounts on Medtrx and DMC.




**GMH-1**  
Glove & Mask box holder. Made from lightweight aluminum. Adjusts to many sizes of boxes. Mounts to railing on cart.



**CW-1**  
Cord wrap. Usually used with outlet strip (OS-1) sold separately.




**LT-1**  
Flexible magnetic light for use on steel carts. Requires two AA batteries. Uses highly efficient LED bulb (included).




**CHT-1**  
Chart Holder. Mounts to accessory railing. 13"W x 4"D x 8.5"H



**UR-1, UR-2, UR-3**  
Utility receptacles, available in three sizes, UR-1 (0.85 liters), UR-2 (1.5 liters), UR-3 (2.8 liters). Requires accessory rail.



**FT-1**  
Fluid/equipment storage tray system for Medtrx/Halftr carts. Includes two trays. Aluminum carts purchase FTA-1. Each tray measures 17.52"W x 4.04"D x 4.5"H



**TRIM-MT**  
Custom color your cart to match you décor. We offer four colors of drawer trim to choose, Blue, Green, Pink, or Violet.



**UR-4L**  
Utility receptacle, 4" opening, 11" Tall



**WNH-1**  
Wire notebook holder. Mounts to accessory railing.

Please use the following part numbers. The last letter will specify color

Half-sized Drawers: TRIM-MTH(B)



Full-sized Drawers: TRIM-MTF(B)

B=blue, G=green, P=pink, V=Violet

TRIM-MTFB, TRIM-MTFG, TRIM-MTFP  
TRIM-MTFV, TRIM-MTHB, TRIM-MTHG  
TRIM-MTHP, TRIM-MTHV

# Medtrx/Halftrx Cart Configurations

Set-up MED-TRX 1.5 (full size) and HALF-TRX 1.0 (compact) medication carts. The chart below indicates the required spaces for each item. (Side A is crossed out on the HALF-TRX 1.0 configuration).



Half size drawers.

Full size drawers and cassettes.

Med-Trax Item #	Description	Spaces Required
DRW-3H	3" Half Size Drawer	2
DRW-4H	4" Half Size Drawer	3
DRW-7H	7" Half Size Drawer	5
DRW-9H*	9" Half Size Drawer	6 2/3
DRW-10H	10" Half Size Drawer	7
DRW-3F	3" Full Size Drawer	2
DRW-4F	4" Full Size Drawer	3
DRW-6F	6" Full Size Drawer	4
DRW-7F	7" Full Size Drawer	5
DRW-9F*	9" Full Size Drawer	6 2/3
DRW-10F	10" Full Size Drawer	7
CASS-3MTC	2 Tier Cassette w/3" bins	4
CASS-5MTC	2 Tier Cassette w/5" bins	4
CASS-7MTC	2 Tier Cassette w/7" bins	4

A completed set-up sheet is required with orders on custom MED-TRX and HALF-TRX Carts. Please contact Customer Service or your Sales Representative if you require assistance with set-up. Custom carts are not returnable.

### STORING NARCOTICS?

Individual locks can be added to medication cart drawers (not available for 3"H). Drawers are keyed separately from cart to add extra security. Ask your sales representative or call our customer service department at 1-800-833-4419 for more information.

Full Drawers: 16"D x 17.25"W  
 Half Drawers: 16"D x 9.125"W

\*Must be ordered in sets of three.

**MEDTRX 1.5**  
 Dimensions (includes casters)  
 Overall shell size: 39"W x 23"D x 46"H

**Halftrx 1.0**  
 Dimensions (includes casters)  
 Overall shell size: 27"W x 23"D x 46"H



NL-MT Narcotics Lock