



WATERLOO
HEALTHCARE
Better Carts for Better Outcomes

CUSTOM CART CONFIGURATIONS



Waterloo Healthcare recognizes that there are unique and specific needs in most healthcare settings. Our product catalog showcases the most popular cart configurations, however we don't want to limit your options.

Design your own configuration with our easy to follow charts!

ELITE ALUMINUM & TRADITIONAL STEEL

Cart Configurations



STANDARD WIDTH

29" Width x 24.5" Depth
(includes stabilizer base)

Drawer Interior Dimensions
22" Wide x 16.5" Deep

JUNIOR WIDTH

23" Width x 24.5" Depth
(includes stabilizer base)

Drawer Interior Dimensions
16" Wide x 16.5" Deep

*Don't see the cart dimensions you need shown above?
Contact Customer Service for more options!*

BUILDING YOUR PART NUMBER IS AS EASY AS 1, 2, 3!

SELECT CART OPTIONS

- Choose your cart width and height:
Width: STANDARD, JUNIOR OR WIDE
Height: SHORT, MEDIUM, TALL OR X-TALL
- Select any combination of 3", 6", 9" OR 12" drawers that fit within your selected carts configurable drawer space
- Select drawer front colors

Each Alpha or Numeric Digit has a meaning
EXAMPLE:

UTGKA-33669-RED



SECTION 1

1st Letter (Cart Width)	2nd Letter (Cart Height)	3rd Letter (Exterior Color)	4th Letter (Lock System)	5th Letter (Cart Material)
U = Standard	S = Short	B = Elec. Blue	K = Key	A = Aluminum
J = Junior	M = Medium	G = Gray	L = Lever	U = Steel
W = Wide	T = Tall	R = Red	S = Simplex	
	X = X-Tall		E = Electronic	
			P = Proximity	

SECTION 2

1st Number (Height of 1st Drawer)	2nd Number (Height of 2nd Drawer)	3rd Number (Height of 3rd Drawer)	4th Number, etc. (Height of 4th Drawer)
--------------------------------------	--------------------------------------	--------------------------------------	--

Totals for Section 2 must equal total Configurable Drawer Space for your Cart Size

Short Cart = 21" Total
 Medium Cart = 24" Total
 Tall Cart = 27" Total
 X-Tall Cart = 30" Total

SECTION 3

3 letters signify color choice for drawer fronts
 (Colors to choose from are listed below)

