



COLD STORAGE VACCINE CART



UXGKA-36H-DKB (Cart Only)
79-2105 (Fridge Only)



- ACCUCOLD 1 cu. ft. refrigerator
- Intelligent microprocessor digital temperature controller
- Adjustable operating control range from +2°C to +10°C
- Digital display of the min/max temperature (viewable in Celsius or Fahrenheit)
- Audible/visual alarms for high/low temperature excursions
- Key locking front door and hospital grade plug come standard
- 115 VAC/60 HZ Only

Vital supplies are always close at hand with our Cold Storage Vaccine Cart. The built in refrigerator allows you to store vaccines at recommended dosing temperatures!

Two key locking drawers allow for additional storage, and two pull out side shelves offer extra workspace. Aluminum construction makes the cart lightweight and easy to maneuver.

DIMENSIONS

CART: 29"(w) x 24.5"(d) x 45"(h)

REFRIGERATOR:

17.5"(w) x 19.75"(d) x 21.5"(h)

RECOMMENDED ACCESSORIES



Sharps Bracket
SB-2L



Accessory Panel
AMP-1



Waste Basket
BSK-2



Horizontal Rail System
BRS-2



Tilt Bins
TB-5



Glove Box Holder
BH-3