



Waterloo
Healthcare Medical Carts

MALIGNANT HYPERTHERMIA CART



UXGKA-36H-DKB (Cart Only)
79-2105 (Fridge Only)



Vital supplies are always close at hand with our Malignant Hyperthermia cart. The built in refrigerator holds up to six 1000ml saline bags and still has room for other items!

Two key locking drawers allow for additional storage, and two pull out side shelves offer extra workspace.

Aluminum construction makes the cart lightweight and easy to maneuver.

DIMENSIONS

CART: 29" (w) x 24.5" (d) x 45" (h)

REFRIGERATOR: 17.5" (w) x 19.75" (d) x 21.5" (h)

- ACCUCOLD 1 cu. ft. refrigerator
- Intelligent microprocessor digital temperature controller
- Adjustable operating control range from +2°C to +10°C
- Digital display of the min/max temperature (viewable in Celsius or Fahrenheit)
- Audible/visual alarms for high/low temperature excursions
- Key locking front door and hospital grade plug come standard



Waterloo
Healthcare Medical Carts

800.833.4419
waterloohealthcare.com