



WATERLOO HEALTHCARE

Better Carts for Better Outcomes

ALL-IN-ONE PHLEBOTOMY CARTS

Our all-in-one phlebotomy cart is a secure mobile solution for storing and transporting phlebotomy supplies throughout the hospital, and built to minimize trips back to the lab for fresh supplies. This cart offers a variety of storage compartments and tilt bins. With it's unique design, customers have also utilized this cart as an IV cart, Suture cart and more.

PAPER HOLDER

Allows you to organize specimen slips and other documents

UNIVERSAL GLOVE BOX

Adjustable to accommodate a single glove box.

REMOVABLE PLASTIC TOP

Facilitates cleaning and allows for easy replacement if necessary

KEY LOCK AND GATE LOCK

Secures 4 doors to the left of the tilt bins, keeping contents locked

INTERNAL LOCKABLE SHARPS COMPARTMENT

Allows for convenient and safe disposal of needles

Accommodates a 12 QT Kendall 85221R sharps container (not included)

DRAWERS WITH REMOVABLE DIVIDER TRAYS

Trays can be moved to facilitate your procedure and provide quick exchange when replenishing supplies

TILT OUT TRASH BIN

Provides convenient waste disposal

TILT BINS

Keep supplies visible and allow for rapid accessibility

STEERING CASTER

Facilitates, maneuverability and two locking casters help keep the Cart in place

MTWA-34696-WHT
OVERALL DIMS:
27"W x 24.5"D x 45"H

MSWA-3469-WHT
OVERALL DIMS:
27"W x 24.5"D x 37.25"H

**AVAILABLE IN TWO CONFIGURATIONS
(TALL & SHORT)**