



Waterloo
Healthcare Medical Carts

QUICK SHIP PRODUCTS

Waterloo Healthcare understands there are times when our customers require products sooner than our standard lead time. When that's the case, we offer a selection of common configurations ready for quick shipment.

These carts are usually ready to ship within a week or less!

**Subject to availability*



USRLU-3369-RED

STEEL UNICART 4 DRAWER ALL RED,
LEVER LOCK FOR TAMPER EVIDENT
PLASTIC SEALS

DRAWERS CONFIGURATION:

- (2) Small: 3" H
- (1) Medium: 6" H
- (1) Large: 9" H

TOTAL DRAWER SPACE: 21"
DRAWERS DIMS: 22"W x 16.5"D
OVERALL DIMS: 29"W x 24.5"D x 36"H



UTRLU-333369-RED

STEEL UNICART 6 DRAWER ALL RED,
LEVER LOCK FOR TAMPER EVIDENT
PLASTIC SEALS

DRAWERS CONFIGURATION:

- (4) Small: 3" H
- (1) Medium: 6" H
- (1) Large: 9" H

TOTAL DRAWER SPACE: 27"
DRAWERS DIMS: 22"W x 16.5"D
OVERALL DIMS: 29"W x 24.5"D x 42"H



Waterloo
Healthcare Medical Carts

800.833.4419
waterloohealthcare.com



ESGKU-399-LTG

ECONOMY CART 3 DRAWER ALL LIGHT GRAY, KEY LOCK
DRAWERS CONFIGURATION:
(1) Small: 3"H
(2) Large: 9"H

TOTAL DRAWER SPACE: 21"
DRAWERS DIMS: 16.25"W x 16.5"D
OVERALL DIMS: 23"W x 24.5"D x 36"H



ESGKU-3369-LTG

ECONOMY CART 4 DRAWER ALL LIGHT GRAY, KEY LOCK
DRAWERS CONFIGURATION:
(2) Small: 3"H
(1) Medium: 6"H
(1) Large: 9"H

TOTAL DRAWER SPACE: 21"
DRAWERS DIMS: 16.25"W x 16.5"D
OVERALL DIMS: 23"W x 24.5"D x 36"H



ESGKU-33366-LTG

ECONOMY CART 5 DRAWER ALL LIGHT GRAY, KEY LOCK
DRAWERS CONFIGURATION:
(3) Small: 3"H
(2) Medium: 6"H

TOTAL DRAWER SPACE: 21"
DRAWERS DIMS: 16.25"W x 16.5"D
OVERALL DIMS: 23"W x 24.5"D x 36"H



UTGKU-3699-YEL

STEEL UNICART 4 DRAWER LIGHT GRAY with YELLOW DRAWERS, KEY LOCK
DRAWERS CONFIGURATION:
(1) Small: 3"H
(1) Medium: 6"H
(2) Large: 9"H

TOTAL DRAWER SPACE: 27"
DRAWERS DIMS: 22"W x 16.5"D
OVERALL DIMS: 29"W x 24.5"D x 42"H



MSWA-3469-WHT

PHLEBOTOMY CART 4 DRAWER ALL WHITE, GATELOCK
HALF-SIZE DRAWERS CONFIGURATION:
(1) Small: 3"H
(1) Small: 4"H
(1) Medium: 6"H
(1) Large: 9"H

TOTAL DRAWER SPACE: 22"
OVERALL DIMS: 27"W x 24.5"D x 37.25"H



MTWA-34696-WHT

PHLEBOTOMY CART 5 DRAWER ALL WHITE, GATELOCK
HALF-SIZE DRAWERS CONFIGURATION:
(1) Small: 3"H
(1) Small: 4"H
(1) Medium: 6"H
(1) Large: 9"H

FULL-SIZE DRAWER CONFIGURATION:
(1) Medium: 6"H

TOTAL DRAWER SPACE: 28"
OVERALL DIMS: 27"W x 24.5"D x 45"H